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# neo.lign A, neo.lign P

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## Instructions for use

# neo.lign A, neo.lign P

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

### 1. Indication range

neo.lign A (Anterior) and neo.lign P (Posterior) are artificial teeth made of acrylic material for the fabrication of dental restorations in the dental laboratory. The shapes (designs) and shades of the neo.lign denture teeth are identical with those of the novo.lign veneers and are particularly suited for the fabrication of combined restorations, implant-supported dentures and for crown and bridgework. The individual, prefabricated tooth design allows for fast fabrication of dentures and veneers and even enables customization.

### 2. Processing

The teeth to be replaced or to be veneered are selected from the range of various shades and design whilst taking the physiognomic characteristics of the patient or esthetic planning into account. The shade is selected with a shade guide, such as the classical A-D shades according to VITA. A design chart is available for selecting the design and the tooth size. The information on the shade and the design is printed on the tooth bar.

The prefabricated teeth are adapted to the specific patient situation by milling (using tungsten carbide burs). The teeth are set up in a wax base; based on the esthetic try-in, dentist and patient can decide on major aspects of esthetics and function and the tooth shade, tooth design and set-up of the teeth. Depending on the systems in use, the wax set-up is invested in plaster (for example) and the wax is boiled out. neo.lign A and neo.lign P denture teeth form a chemical bond and bond to denture base materials based on heat- or cold-curing polymethyl methacrylates (PMMA-based) and thermoplastic, PMMA-based denture base materials. To achieve perfect bonding, the basal surfaces of the teeth are roughened, for example using a tungsten carbide bur, or removed from the flask and sandblasted with 110 µm aluminium oxide at a pressure of 2-3 bars. The bond strength can be increased by applying a bonding agent, for example Haftconnector bredent REF 520 0029 2, to the sandblasted or roughened surfaces. The use of a bonding agent, however, is not required unless the use is prescribed in the instructions for use of the respective denture base material. Then the denture base resin is pressed or injected in accordance with the instructions for use.

Subsequently, crea.lign can be used both for customizing the teeth and the denture base; please observe the information in the instructions for use. After finishing, the complete denture including the teeth is polished.

### 3. Side effects and risks

Limited cases of allergic reactions may result from grinding the teeth owing to the release of PMMA. Please notify bredent immediately, if other side effects and risks occur.

### 4. Cleaning of dentures

Do not use reactive denture cleaning agents for cleaning dentures. Use standard tooth pastes and denture cleaning brushes for cleaning.

### 5. Symbols

Lot number

### 6. Additional information

The information contained in these instructions for use is always updated according to the latest knowledge and experience. Therefore we recommend to read the instructions for use again before using a new package. These non-binding user recommendations are based on our own experience. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since bredent does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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